

〈For Research Use Only〉 Code No. 296-63801 (1,300 Tests)

Wako LabAssay™ Phospholipid (Choline Oxidase · DAOS method)

〔Introduction〕

Phospholipids are a class of lipids, and a major component of all biological membranes, along with glycolipids, cholesterol and proteins.

LabAssay™ Phospholipid is based on an enzymatic methods using *N*-ethyl-*N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS) as a blue pigment. This kit is used for the quantitative determination of phospholipids in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

〔Kit contents〕

(1)	Buffer (50mmol/L Good's buffer, pH7.5)	50mL	8vials
(2)	Chromogen Substrate (when reconstituted)	For 50mL	8vials
	Phospholipase D	0.47units/mL	
	Choline Oxidase	2.0units/mL	
	Peroxidase	4.2units/mL	
	<i>N</i> -Ethyl- <i>N</i> -(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS)	0.77mmol/L	
	4-Aminoantipyrine	0.24mmol/L	
	Ascorbate Oxidase	3.9units/mL	
(3)	Standard Solution (Choline Chloride 54mg/dL (Corresponding to 300mg/dL Phospholipid))	10mL	2vials

〔Materials and apparatuses to be prepared〕

- 96-well microplate (transparent type)
- Micropipette
- Incubator maintaining at 37°C*
- Plate mixer*
- Microplate reader with a 600nm wavelength filter
(* if the microplate reader is not equipped.)

〔For test tube method〕



- Test tube
- Pipette
- Spectrophotometer or colorimeter with 600 nm wavelength filter

〔Assay principle〕


Phospholipids (lecithin, sphingomyelin and lysolecithin) in a sample are hydrolyzed to choline in a reaction catalyzed by phospholipase D. Choline so formed is oxidized by choline oxidase in a reaction that produces hydrogen peroxide. The hydrogen peroxide produced causes DAOS and 4-Aminoantipyrine to undergo a quantitative oxidative condensation catalyzed by peroxidase (POD), producing a blue pigment. The amount of

-1/8-


phospholipids contained in the sample is determined by measuring the absorbance of the blue color.



Choline Oxidase



POD



〔Preparation of reagents to be used〕

① Color Reagent :

Prepare Color Reagent by Dissolving 1 vial of Chromogen substrate (for 50mL) to 50mL of Buffer. After reconstitution, the solution should be stored at 2-10°C and used within one week.

② Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard.

No.	Phospholipid Standard	Distilled or deionized water	Sample volume	Phospholipid Concentration
1	100 μL	100 μL	2 μL	150 mg/dL
2	Undiluted	—	2 μL	300 mg/dL
3	Undiluted	—	4 μL	596.1 mg/dL* ¹

*1 The test sample volume is usually 2 μL, but 4 μL is taken in this case. The phospholipids concentration must be corrected accordingly as indicated in the table above.

〔Procedure〕

(1) Assay in a microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
Sample	Serum 2 μL	Standard solution 2 μL	—
Color reagent	300 μL	300 μL	300 μL
	Mix well and incubate at 37°C for 5 min. Measure the absorbance at 600nm* ² of the test sample and standard solution with the blank solution as the control.		

*2 In two wavelength assay, measure using Main wavelength 600nm/Sub wavelength 700nm.

-2/8-


(2) Assay in a test tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Sample	Serum 0.02 mL	Standard solution 0.02 mL	—
Color reagent	3 mL	3 mL	3 mL
	Mix well and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600nm Spectrophotometer : 600nm*2		

*2 In two wavelength assay, measure using Main wavelength 600nm/Sub wavelength 700nm.

[Standard curve] [assay in a microplate]



Microplate reader : SAFIRE (TECAN)

[Performance]

(1) Sensitivity [assay in a test tube]

- The absorbance is below 0.14, when measuring purified water as a sample.
- The absorbance is 0.12 to 0.58, when measuring 300mg/dL phospholipid as a sample.

(2) Specificity

- The phospholipid concentration is less than $\pm 10\%$, when measuring the known concentration of control serum as a sample.

[Usage Notes]

(1) Sample

- Ascorbic acid, bilirubin and hemolysis may not significantly affect the assay.
- This method measures lecithin, sphingomyelin, lysolecithin but not cephalin.

(2) Notes on the assay

- Do not use the reagents described above in any procedures other

than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.

- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- The reaction completes in about 3 min. It is enough to incubate at 37 °C for 5 min. When the incubation is continued, the absorbance gradually decreases. Therefore, the incubation should be done within 15 min.
- Color level of the reaction solution should show very little change within 1 hour at room temperature.
- This kit is for research use only. Not for diagnostic use.

(3) Safety precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette should be used with a safety pipette filler.

(4) Waste

- The waste should be processed appropriately according to the law.
- All the devices including reagents and vials that come into contact with specimens should be considered potentially infectious.

Expiration date : 12 months after the manufacture

Storage : Store at 2 ~ 10°C

Package : 1,300 tests

[Reference]

1. Takayama, M., Itoh, S., Nagasaki, T. and Tanimizu, I. : *Clin. Chim. Acta.*, **79**, 93 (1977).

Wako Pure Chemical Industries, Ltd.

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.wako-chem.co.jp>

Wako Chemicals USA, Inc.

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : +1-804-271-7677
Facsimile : +1-804-271-7791
<http://www.wakousa.com>

Wako Chemicals GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-3111-0
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>